

1 Kb YEA Ladder DNA Marker II

Ready-to-use DNA ladder

Recommended storage condition: -20 °C

Loading volume: 5~10 μ l/well

Size range: 0.25~10 kb

Concentration: 60 ng/ μ l

Package size: 60 μ g/1 mL

1 Kb DNA Ladder, ready-to-use

The Yeastern Biotech's 1 Kb YEA DNA Ladder Marker II (FYD008-1ML) contains 13 discrete fragments ranging from 0.25 kb to 10 kb. It is supplied in a ready-to-use format containing blue tracking dyes, and mixed at varied concentrations.

For best results, load 10 μ L of the marker per well. It allows sizing and concentration estimation of DNA fragments on agarose gels generated by PCR or restriction digestion.



www.yeastern.com


yeastern@yeastern.com

t / 886-2-26953922 f / 886-2-26953979

Cat. No. FYD008-1ML

1 Kb YEA Ladder DNA Marker II

Ready-to-use DNA ladder



Loading volume: 10 μ l/well
1.0%, 1 \times TAE agarose gel



Yeastern Biotech Co., Ltd.